

Freedom Excel Superlite Self-Propelled Wheelchair - 20"

SKU: MWS452180

\$697



Features and Benefits:

- Folding backrest for easy transportation
- Flip up removable padded armrest
- Park brake below seat level when engaged for easy side transfer
- Outwards swinging removable legrest
- Height Adjustable Plastic Footplates
- Airless "soft ride" polyurethane tyres
- Quick release removable spoke wheels
- Double-braced aluminium construction
- Spill proof nylon fabric upholstery
- Footrest Length: 300-400mm
- Folded Width: 290 mm

Specifications:

Code	MWS452180	Safe Working Load	110kg
Description	Freedom Excel Superlite Self-Propelled Wheelchair – 18" (46cm)	Unit Weight	12.7kg
Total Length	1070mm	Seat Depth	430mm
Total Width	690mm	Seat Width	510mm
Total Height	950mm	Seat to Floor Height	500mm
Seat to Footrest	500mm	Backrest Height	450mm