

# Freedom Excel Superlite Transit Wheelchair - 18"

SKU: MWS452110 (AWC137T)

**\$697**



## Features and Benefits:

- Folding backrest for easy transportation
- Handbrake with park locking mechanism connected to hub brakes
- Double-braced aluminium construction
- Drum hub brake system
- Height adjustable plastic footplates
- Outwards swinging removable legrest
- Airless "soft ride" polyurethane tyres
- Spill proof nylon fabric upholstery
- Footrest Length: 300-400 mm
- Folded Width: 290 mm

## **Specifications:**

<b>Code</b>	MWS452110	<b>Safe Working Load</b>	110kg
<b>Description</b>	Freedom Excel Superlite Transit Wheelchair – 18" (46cm)	<b>Unit Weight</b>	11.8kg
<b>Total Length</b>	920mm	<b>Seat Depth</b>	430mm
<b>Total Width</b>	640mm	<b>Seat Width</b>	460mm
<b>Total Height</b>	950mm	<b>Seat to Floor Height</b>	500mm
<b>Seat to Footrest</b>	500mm	<b>Backrest Height</b>	450mm