

Freedom Bath Transfer Bench – HBA424

SKU: BTS153100

\$263



Features and Benefits:

- Removable backrest can be attached to either side of the bench
- Detachable armrest
- Height adjustable legs. Adjustment pins are spaced at 25mm increments on three extension legs and 12.5mm increments on the fourth to maximise stability
- Thicker, more durable powder coated aluminium frame to prevent scuffs and scratches
- Heavy duty sure grip suction rubber tips
- Purchased separately: Extension Legs (GSP650100)

Specifications:

Code	BTS153100	Seat Width	680mm
Description	Freedom Bath Transfer Bench – HBA424	Seat Depth	410mm
Overall Width	785mm	Seat Height	490mm – 640mm
Overall Depth	520mm	Safe Working Load	200kg
Overall Height	825mm – 975mm	Warranty	1 Year