

Aspire Vogue Carbon Fibre Seat Walker

SKU: WAF705350

\$750



Features and Benefits

- The stylish European designed Aspire Vogue Walkers are a lightweight, versatile range with a sturdy frame and easy to fold locking system, suitable for compact transportation when folded
- Smooth height adjustment with AUTOLOCATE - Graduated and numbered upright shaft automatically locates to preset heights, making adjustment simple and accurate
- Super Lightweight 5 kg Carbon Fibre Frame
- Luxury Look and Feel - European design, profiled tubing, moulded fittings and chic colour palette
- Lockable Handbrakes - For increased user safety when walking and when sitting
- Lightweight and Ergonomic - Easy to push with a small turning circle and low rolling resistance
- Unbreakable TPE Super Soft Wheels - High durability with a smooth ride for increased comfort
- Locks when Folded - Makes lifting the frame easier and minimises storage space required
- Removable Utility Bag - Stylish and simple to remove for shopping

Specifications

Code	WAF705350	Seat Depth	190mm
Description	Vogue Walker – Carbon Fibre	Seat Width	460mm
Overall Height	695mm – 920mm	Seat Height	550mm
Overall Length	680mm	Wheel Size	8"
Overall Width	600mm	Total Weight	5kg
Folded Width	225mm	Bag SWL	5kg