



BME

BRIBIE MEDICAL EQUIPMENT

*Local Family
Operated Business*

PRODUCT SPECIFICATIONS SHEET

SKU: WAF705340 (Short)

SKU: WAF705350 (Medium)

SKU: WAF705370 (Tall)



Sapphire
Short



Black
Medium / HD



Emerald Green
Tall

Aspire Vogue Carbon Fibre Seat Walker

PRODUCT OVERVIEW

- Lightweight carbon fibre frame
- Ergonomic design with excellent manoeuvrability
- Foldable frame with locking mechanism for transport and storage
- Locking handbrakes for safety
- Soft, durable wheels for a smooth ride
- Removable storage bag
- Includes cane holder



SPECIFICATIONS

	Short (Sapphire)	Medium (Black)	Tall (Emerald)
Total Length	655 mm	655 mm	655 mm
Total Width	600 mm	615 mm	600 mm
Total Height (Adjustable)	715 - 830 mm	725 - 940 mm	780 - 1030 mm
Folded Width	230 mm	230 mm	230 mm
Seat Depth	240 mm	240 mm	240 mm
Seat Width	460 mm	460 mm	460 mm
Seat Height from Ground	500 mm	550 mm	620 mm
Wheel Size	8"	8"	8"
Unit Weight	5.9 kg	6.2 kg	6.35 kg
Bag Safe Working Load	5 kg	5 kg	5 kg
Safe Working Load	150 kg	150 kg	150 kg