

Aspire Vogue Carbon Fibre Seat Walker - TALL

SKU: WAF705360 (Candy Red)

\$760



Features and Benefits

- **Super lightweight** - Carbon fibre frame.
- Now supplied with bonus cane holder (ref. image)
- **Luxury look and feel** – Danish design, profiled tubing, moulded fittings, ergonomic handles and refined colours.
- **Lightweight and ergonomic design** with small turning circle for excellent manoeuvrability.
- **“AutoLocate”** - Smooth height adjustment that automatically locates to pre-set heights, making adjustment simple and accurate.
- **Locks when folded** – makes transportation and storage easier.
- **Advanced locking handbrakes** – for increased user safety when walking and sitting.
- **Unbreakable TPE super soft wheels** – offers high durability and a smooth ride for increased comfort.
- **Removable utility bag** – simple and stylish to remove for shopping. Can hold up to 5 kg.

Specifications:

Code	WAF705360	Seat Width	460mm
Description	Aspire Vogue Carbon Fibre Seat Waker - TALL	Seat to Floor Height	620mm
Total Length	680mm	Colour	Candy Red
Total Width	605mm	Safe Working Load	150kg
Total Height	770mm – 1040mm	Unit Weight	6kg