

Omega PA1 Self-Propelled Children's Wheelchair

SKU: 62005

\$796



Features and Benefits:

Small Adjustable Self-Propelling Wheelchair designed for children or smaller users.

- Adjustable, removable arm rests.
- Swing in/out leg rests.
- 22" quick-release rear wheels.
- Rear anti-tip wheels for added security.
- User weight capacity 125kg.

Specifications:

Approximate Total Product Weight	13.7kg	Overall Width when Opened	565mm
Overall Width when Closed	300mm	Approximate overall depth	600mm
Seat Width	355mm	Seat Depth	355mm
Seat Height (to floor)	445mm-470mm	Front Castor	20cm x 2.5cm