

Pathrider 140XL

SKU: PR140XL-4Y (Yellow)

SKU: PR140XL-4S (Silver)

SKU: PR140XL-4R (Red)

\$6,497



International seat shown

Features and Benefits:

- Top of the range sleek and sporty powered mobility scooter
- Large 13" pneumatic tyres to handle rugged outdoor terrain.
- Luxurious suspension seat with arm rests, head rest and lumbar support for stability
- Front and rear suspension and sealed hydraulic brakes.
- Automatic full lighting system including directional adjustment. Long lasting high intensity LED kerb light, headlight, taillight, and indicators
- Back lit battery gauge
- Wrap around easy drive tiller.
- Tiller mounted battery charger port for off board charging.
- Unit Weight: 118 kg

Warranty:

Two-Year Limited Warranty:

- Structural Frame Components (platform, fork, seat post, frame welds)
- Drivetrain (Differential, motor, brake)

One-Year Limited Warranty:

- Your Pathrider 140XL is fully guaranteed for 12 months from the date of purchase against faults arising due to defects in manufacture or materials. This warranty does not detract from, but is in addition to your legal rights.
- All electric parts (controllers, battery chargers). Any attempt to open or dismantle these items renders the guarantee void on this item.

Not Covered Under Warranty:

- Normal wear and tear (tyres, belts, bulbs, upholstery, plastic shrouds, motor brushes, fuses, batteries)
- Damage to the product caused by misuse or accident


BATTERIES:

- 6-month Warranty
- Gradual deterioration in performance due to being left in a discharged state, left in cold conditions for long periods of time, or worn out through heavy use is not covered.

197 First Ave, Bongaree, QLD, 4507 – bribie@medicalequipmenthire.com.au – (07) 3408 0398



Specifications:

| | |
|---|---|
| Weight Capacity | 181 kgs (400 lbs.) |
| Maximum Speed¹ | Up to 10 km/h (6.2 mph) |
| Ground Clearance² | 13.33 cm (5.25") |
| Turning Radius² | Right: 210.8 cm (83") Left: 208.3 cm (82") |
| Overall Length² | 142.8 cm (56.25") |
| Overall Width² | 69.2 cm (27.25") |
| Front Tyres | 10 x 33 cm (4" x 13") pneumatic |
| Rear Tyres | 10 x 33 cm (4" x 13") pneumatic |
| Range Per Charge^{2,4} | Up to 35 km (21.70 miles) (with 75 AH batteries) |
| Total Weight Without Batteries | 117.70 kg (259.5 lbs.) |
| Heaviest Piece When Disassembled | 93.21 kg (205.5 lbs.) (mid-frame) |
| Standard Seating | Type: Pillow-top Weight: 21.77 kg (48 lbs.) Material: Black vinyl Dimensions: Width: 49.5 cm (19.5") Depth: 44.5 cm (17.5") |
| Motor Output | 24-volt DC motor, rear-wheel drive, sealed transaxle |
| Brakes | Regenerative and electromechanical |
| Battery Requirements^{3,5} | (2) 12V deep-cycle Size: 70 AH (standard) Weight: 70 AH: 25 kg (55 lbs.) each |
| Battery Charger | Off-board, 8 amp |
| Suspension | Front and rear |
| Seat-To-Floor Height | 72.39 - 68.58 cm (28.5" - 27") |
| Colour Options |  |