



BME

BRIBIE MEDICAL EQUIPMENT

Local Family

Operated Business

PRODUCT SHEET

SKU: VR205SP

Ashford Toilet Surround



PRODUCT OVERVIEW

- Toilet frame surround for help getting on and off the toilet
- Moulded arms for improved grip
- Height and width adjustable
- Double front bar for strength

SPECIFICATIONS

Overall Height	660 – 810 mm	Overall Depth	450 mm
Overall Width	580 mm	Unit Weight	4.3 kg
Width Between Arms	475 – 525 mm	Max. Weight Capacity	190 kg
Material	Epoxy-coated Steel		