



BME

BRIBIE MEDICAL EQUIPMENT

*Local Family
Operated Business*

PRODUCT SHEET

SKU: BTT155100

Aspire Aluminium Toilet Surround

PRODUCT OVERVIEW

- Raised armrests for safe and stable transfers
- Push-button adjustable legs for customizable height
- Independent leg adjustment for uneven surfaces
- Non-slip rubber feet for safety on wet floors
- Easy to clean and maintain
- Suitable for use with standard toilets

SPECIFICATIONS

Material	Aluminium	Width Between Arms	450 mm
Overall Height	535 - 660 mm	Unit Weight	3 kg
Overall Width	500 mm	Safe Working Load	175 kg
Overall Depth	435 mm		