

# Aspire Vogue Lightweight 2 Seat Walker

SKU: WAF705310CH (Champagne)

SKU: WAF705310RE (Red)

SKU: WAF705310LE (Leopard)

SKU: WAF705310FL (Flower)

SKU: WAF705310ZB (Zebra)

**\$497**

**\$589**



## **FEATURES AND BENEFITS:**

This modern seat walker has a clean and crisp frame design, making it a stand out when going for a walk or browsing the shops.

- Lightweight aluminium body improves manoeuvrability and portability
- Ergonomic Design Offers personalised height adjustment which promotes correct posture
- Compact turning Circle for improved mobility in tight spaces
- Simple folding mechanism locks for easy transportation and storage
- Minimised storage space with folded width 220 mm
- Handbrakes Lock for user safety and stability
- Removable bag carries up to 5kg of essentials
- Deep seat and wide backrest for improved comfort and support
- Bonus Walking cane holder

Total Length	Total Width	Total height	Seat Width	Safe Working Load	Unit Weight	Seat to Floor Height	Frame Weight
655mm	620mm	685-920mm	465mm	150kg	6.45kg	540mm	5.3kg

