

Alivio Donatello Dual Motor Lift Recliner – Petite

SKU: KA555

\$3,375



Features and Benefits:

- Available in Merino Steel (3 roll backrest and lateral backrest) and Canyon (lateral backrest only)
- Two stage lift function
- Lay-flat chair.
- Footrest extension
- Easy-to-use hand controller
- Ultra-durable stylish fabric
- Lumbar support to accommodate sitting posture.
- Power headrest to lift or lower head position.

Specifications:

Overall Depth	978mm	Seat Height	457mm
Overall Width	851mm	Reclined Length	1830mm
Overall Height	1041mm	Product Weight	70kg
Seat Depth	457mm	Max User Weight	158kg
Seat Width	457mm	Colour	<ul style="list-style-type: none"> - Canyon (Lateral Backrest only) - Merino Steel (Standard 3 roll and lateral backrest)