

Freedom Bariatric Walking Frame

SKU: WAF752162 (Black)

\$142



Features and Benefits:

- PVC handgrips
- Strong steel construction
- Replaceable skid resistant rubber tips
- Height adjustable extension legs
- Easy-fold mechanism for storage and transportation
- Folded width: 95 mm

Specifications:

Code	WAF752162 (Black), RBE804		
Description	Freedom Bariatric Side-Folding Walking Frame	Safe Working Load	270kg
Total Length	460mm	Unit Weight	3kg
Total Width	740mm		
Total Height	760-1010mm		