



BME

BRIE MEDICAL EQUIPMENT

*Local Family
Operated Business*

PRODUCT SHEET

SKU: BTT163200

Freedom Over Toilet Aid

PRODUCT OVERVIEW

- Non-flexible seat and lid easily removed for cleaning
- Anti-slip suction rubber tips
- Height adjustable legs
- Powder coated aluminium frame
- Splash guard included
- Old Product Code: HBA390A



SPECIFICATIONS

Total Length	510 mm	Seat to Top of Armrest	190 mm
Total Width	550 mm	Seat to Floor Height	410 - 530 mm
Seat Depth	370 mm	Unit Weight	4 kg
Seat Width	470 mm	Safe Working Load	130 kg
Material	Aluminium		