

# Freedom Over Toilet Aid – HBA490A

SKU: BTT163200

**\$109**



## Features and Benefits:

- Non-flexible seat and lid easily removed for cleaning.
- Anti-slip suction rubber tips
- Height adjustable legs
- Powder coated aluminium frame.
- Splash guard included.

## Specifications:

<b>Code</b>	BTT163200	<b>Overall Height</b>	605mm – 730mm
<b>Description</b>	Freedom Over Toilet Aid	<b>Seat Width</b>	360mm
<b>Overall Width</b>	550mm	<b>Seat Depth</b>	430mm
<b>Overall Depth</b>	510mm	<b>Seat Height</b>	410mm – 530mm
<b>Safe Working Load</b>	130kg		