

Karma Ergo Lite Deluxe Transit Wheelchair – 16”

SKU: KCKM2512F14-16

\$1,179



Features and Benefits:

- Ultra lightweight, folding, attendant propelled or self-propelled style wheelchair
- S-Ergo patented seating system for added comfort
- Push and pull brake system for 14" rear wheel
- Detachable swing away footrests
- Optional fixed or quick release 20" self-propel wheel
- Armrests designed to allow closer seating to desks
- Top backrest section and frame folds for ease of transport
- S-Ergo patented seating system
- Adjusts and stabilises seating pose
- Double cross braces provide strength and durability

Specifications:

Max User Weight	100kg	Seat Depth	420mm
Overall Length	935mm	Floor to Seat Height	480mm
Overall Width	570mm	Folded Width	280mm
Overall Height	900mm	Folded Depth	635mm
Seat Width	410mm	Armrest Height	215mm
Backrest Height	420mm	Product Weight	8.7kg