

Aspire Monet Dual Action Lift Recline Chair

SKU: CHP227540 (Canyon Blue)

\$2,949



Features and Benefits:

- **Dual Action** - Independent control of leg and recline for infinite positioning
- **Independent Adjustable Headrest** improves comfort whilst watching TV or reading a book
- **Integrated Lift and Recline** functions for easy sit-to-stand
- **Extending Legrest** supports more of the lower legs for enhanced comfort
- **Heated Function** within the backrest provides warmth and comfort during the colder months.
- **Zero Gravity** tilt-in-space positioning creates a feeling of weightlessness and promotes pressure redistribution
- **Backlit Hand Control** with integrated USB charging point
- **Lithium Ion Battery Backup** ensures chair can be used in the event of a power outage or if moved between rooms during the day

Specifications:

Code	CHP227540	Seat Depth	535mm
Description	Aspire Monet Dual Action Lift Recline Chair	Seat Width	560mm
Total Width	890mm	Seat to Floor Height	495mm
Total Height	1065mm	Seat to Footrest	495mm
Safe Working Load	159kg	Backrest Height	570mm