

Aspire Vogue Adventure Walker

SKU: WAF705450 (Black)

\$789



Features and Benefits

- **Heavy-duty** frame and wider seat provide enhanced support for users up to 200kg.
- **European Design** offers personalised height adjustment, promoting correct posture.
- **Large Multi-Terrain Wheels** for durability and stability on and off road. (Front Wheels 11", Back Wheels 9")
- **Ergonomic Handbrakes** lock for improved user safety and stability.
- **Compact Turning Circle** facilitates improved mobility in tight spaces.
- **Simple Folding Mechanism** enables easy transportation and compact storage (245 mm folded width, perfect for the car).
- **Detachable Storage Bag** carries up to 5kg of essentials.
- **Bonus** walking cane holder provides increased flexibility for your mobility needs.
- Unit weight includes accessories, bag and backstrap

Specifications

Code	WAF705450	Safe Working Load	200kg
Description	Aspire Vogue Adventure Walker	Unit Weight	10kg
Total Length	770mm	Seat Width	540mm
Total Width	695mm		
Total Height	900mm – 1035mm		