

Aspire 46cm Weight Bearing Platform Shower Commode – Attendant Propelled

Inc. Closed Seat

SKU: BTC066200

\$1,419



Features and Benefits:

- 8° backrest angle ensures patient comfort
- Improved over toilet clearance
- Integrated access trough improves hygiene whilst maintaining toilet bowl clearance height of 465mm
- 4x Tente total locking castors provide brake access from all angles
- Safety arms prevent falls and swing back out of the way for transfers
- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Easy-slide weight bearing platform footrest slides completely away for transfers

Specifications:

Code	BTC066000	Backrest Height	305mm
Overall Height	995mm	Armrest Height	225mm
Overall Width	535mm	Anterior Seat Height	590mm
Seat Width	460mm	Posterior Seat Height	560mm
Safe Working Load	200kg		