

Aspire XL Seat Walker

SKU: WAF750020SL (Silver)

\$276



Specifications:

Code	WAF750020SL	Seat Width	495mm
Overall Height	850-982mm	Seat Height	560mm
Overall Length	697mm	Backrest Height	850mm
Overall Width	732mm	Folded Width	230mm
Seat Depth	340mm	Unit Weight	8kg
Castor Size	8"	SWL	180kg

Features and Benefits:

- Ultra-soft hand brakes with palm-ball lock effortlessly, offering increased comfort for users with arthritis or wrist/hand joint pain
- Handle height adjusts quickly with hand wheels for optimal sizing
- Rustproof, robust aluminium frame provides durability and stability and whilst remaining very light to manoeuvre and transport
- Easy folding mechanism facilitates storage and transportation
- Anatomically designed, latex-free handles allow comfortable grip especially for users with arthritis
- Curved and padded backrest offers comfortable seated support and flips up/removes easily for folding/storage
- Padded seat flips up easily for access to storage bag (supplied) or for folding
- Swivelling front castor wheels increase manoeuvrability 8" diameter solid, puncture-proof wheels provide increased manoeuvrability over various surfaces
- Easy wipe-clean surfaces
- Wide padded seat and increased 180kg safe working load for large users

