

Aspire Evoke 2 HD Self-Propelled Wheelchair – 61cm

SKU: MWS449550

\$1,129



Features and Benefits:

- A lightweight, highly adjustable manual wheelchair
- The contoured padded overlay offers the user increased comfort and is ideal for those who are in their wheelchair for longer periods of time. The mild contouring assists in positioning and pressure distribution, while the fabric allows air movement, increasing breathability. It is easily removed to allow washing when required
- Swing-in / out, height adjustable and removable foot rests
- Strong, Nylon upholstery with a comfortable cushion and padded backrest.
- Depth adjustable seating with rear angle adjustability
- Quick release 24" solid rear tyres

Specifications:

Seat Depth	400mm - 475mm	Seat Width	610mm
Backrest Height	430mm – 490mm	Seat to Floor Height	440mm – 480mm
Colour	Silver	Material	Aluminium
Safe Working Load	180kg	Unit Weight	16.2kg