



BME

BRIE MEDICAL EQUIPMENT

*Local Family
Operated Business*

PRODUCT SHEET

SKU: 80301214

Etac Cloo Toilet Seat Raiser with Armrests

PRODUCT OVERVIEW

- Tilttable seat assists standing and reduces hip and knee pressure
- Detachable, foldable armrests provide support and flexibility
- Tool-free installation on standard toilets
- Removable seat unit for easy cleaning
- Discreet design suits most bathroom interiors

SPECIFICATIONS

Seat Width	400 mm	Seat Height	60 / 100 / 140* mm
Armrest Height from Seat	210 mm	Opening of Toilet	290 (D) x 220 (W) mm
Width Between Arms	490 mm	Overall Width	600 mm
Unit Weight	4.8 kg	Max. Weight Capacity	150 kg

*To obtain 140 mm (5½") seat height, a set of extra spacers is required.