



BME

BRIBIE MEDICAL EQUIPMENT

Local Family

Operated Business

PRODUCT SHEET

SKU: 80301204

Etac Cloo Toilet Seat Raiser

without Armrests

PRODUCT OVERVIEW

- Tiltable seat reduces hip and knee pressure
- Removable seat unit for thorough cleaning
- Fold-up seat and lid for easy toilet access
- Tool-free installation on standard toilets
- Detachable armrests available separately for added support



SPECIFICATIONS

Seat Width	400 mm	Seat Height	60 / 100 / 140* mm
Unit Weight	3.3 kg	Opening of Toilet	290 (D) x 220 (W) mm
Max. Weight Capacity	150 kg		

*To obtain 140 mm (5½") seat height, a set of extra spacers is required.