



BME

BRIBIE MEDICAL EQUIPMENT

Local Family

Operated Business

PRODUCT SHEET

SKU: BTS118015

Aspire Shower Chair - Wide

PRODUCT OVERVIEW

- Reinforced wide frame for enhanced stability
- Contoured seat with drainage holes for comfort and support
- Corrosion-resistant frame construction
- Raised armrests to assist transfers
- Durable non-slip feet for wet area stability
- Easy-clean design with improved drainage

SPECIFICATIONS

Overall Height	735 - 910 mm	Overall Depth	570 mm
Overall Width	570 mm	Seat Depth	400 mm
Seat Width	440 mm	Seat Height (Adj.)	375 - 550 mm
Width between Arms	525 mm	Unit Weight	3.3 kg
Max. Weight Capacity	200 kg	Material	Aluminium