

Configura Comfort Dual Motor Lift Chair

Available in Beige and Black



From **\$5,293**

Beige - \$5,293

Black - \$5,499



Features and Benefits:

- Designed to provide comfortable seating to users with mild to moderate postural requirements, help reduce the formation of pressure ulcers and help to stand or sit users with limited mobility.
 - 'Reduced shear' recline to maximise pelvic stability
 - Pressure reducing fabric on backrest and armrest as standard
 - Castors for easy manoeuvrability
 - Integrated self-supply overlay for easy integration of existing cushions
 - 2 motors

*Available in Small, Medium, and Large

Specifications: Small

Code	SKU: CR5443 (Black) SKU: CR5422 (Beige)	Armrest Height	180mm
Description	Configura Comfort Dual Motor Lift Chair - Small	Rise height to front of chair STANDARD	546mm
Seat Height	406mm	Rise height to front of chair MEDIUM	622mm
Seat Height with Haks	457mm	Chair to Wall/Furniture	457mm
Seat Width	457mm	Safe Working Load	130kg
Seat Depth Position	457mm	Back Height	700mm
Overall Width	711mm	Overall Height	1092mm
Overall Length Upright Position	940mm	Weight of Product	65kg
Overall Length Recline Position	1702mm		

Specifications: Medium

Code	SKU: CR5444 (Black) SKU: CR5424 (Beige)	Armrest Height	180mm
Description	Configura Comfort Dual Motor Lift Chair - Medium	Rise height to front of chair STANDARD	508mm
Seat Height	457mm	Rise height to front of chair MEDIUM	572mm
Seat Height with Haks	508mm	Rise height to front of chair FLAT	711mm
Seat Width	508mm	Chair to Wall/Furniture	483mm
Seat Depth Position	508mm	Safe Working Load	160kg
Overall Width	762mm	Back Height	700mm
Overall Length Upright Position	940mm	Overall Height	1143mm
Overall Length Recline Position	1702mm	Weight of Product	68kg

Specifications: Large

Code	SKU: CR5445 (Black) SKU: CR5426 (Beige)	Armrest Height	180mm
Description	Configura Comfort Dual Motor Lift Chair - Large	Rise height to front of chair STANDARD	508mm
Seat Height	457mm	Rise height to front of chair MEDIUM	572mm
Seat Height with Haks	508mm	Rise height to front of chair FLAT	711mm
Seat Width	558mm	Chair to Wall/Furniture	483mm
Seat Depth Position	508mm	Safe Working Load	160kg
Overall Width	812mm	Back Height	700mm
Overall Length Upright Position	940mm	Overall Height	1143mm
Overall Length Recline Position	1702mm	Weight of Product	71kg