

Freedom Oval Tube Shower Stool

SKU: HBA407



\$184



Product Description

- Solid plastic seat, perforated for drainage
- Height adjustable legs. Adjustment pins are spaced at 25mm increments on three extension legs and 12.5mm increments on the fourth to maximise stability
- S-bend shaped frame, provides maximum space between the armrests
- Splash guard included
- Thicker, more durable powder coated aluminium frame to prevent scuffs and scratches
- Slip resistant suction rubber tips

Overall Height	Overall Width	Overall Length	Seat Depth	Seat Width	Seat Height	Product Weight	SWL
795mm - 945mm	615mm	605mm	400mm	460mm	440mm – 590mm	3.8kg	200kg

