



BME

BRIBIE MEDICAL EQUIPMENT

*Local Family
Operated Business*

PRODUCT SHEET

SKU: BTS150200

Freedom Shower Chair

PRODUCT OVERVIEW

- Rigid seat made from high quality plastic
- Anti-slip rubber tips
- Height adjustable legs
- Powder coated aluminium frame



SPECIFICATIONS

Overall Height	770 - 895 mm	Overall Length	510 mm
Overall Width	550 mm	Seat Depth	400 mm
Seat Width	470 mm	Seat Height (Adj.)	415 - 540 mm
Seat to Armrest Height	190 mm	Unit Weight	3.6 kg
Max. Weight Capacity	130 kg	Material	Aluminium