

Freedom Shower Seat

SKU: HBA391

\$103



Features and Benefits:

- Frame made from steel tubing.
- Ideal for domestic use.
- Powder coated and corrosion resistance treated.
- Height adjustable extension legs.
- Solid plastic seat with drainage perforations.

Specifications:

Code	HBA391	Seat Width	470mm
Description	Freedom Shower Seat	Seat Depth	400mm
Colour	White	Seat to Floor Height	415mm – 540mm
Total Width	550mm	Safe Working Load	130kg
Warranty	1 year	Unit Weight	3.9kg