

Shoprider ALLROUNDER

9LSPB- Blue

9LSPG- Grey

9SLP- Red



\$4,999



Features and Benefits:

- Easy access control
- User friendly tiller / dash
- USB charger
- Accessory / storage pocket
- Stylish fluted rims
- Double nudge bar
- LED indicators
- High back contoured seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- LED Headlight
- Rear lights
- Hazard lights
- Large battery capacity
- Battery power indicator
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width arm rests
- Removable shopping basket

Specifications:

Overall: 1300mm (L) x 640mm (W) x 1200mm (H)

Turning Circle: 1440mm

Speed: 10km/h

Maximum range: Up to 32km

Safe climb angle: 9°

Maximum climb angle: 12°

Motor: 1.5Hp high torque motor (4 pole)

Batteries: 2 x 12V 50 Amp minimum (supplied)

Charger: 8 Amp

Tyres (Pneumatic): Ø330 x 100mm (4.00-5) – Front + Rear

Product Weight: 99kg – with batteries

Maximum User Weight: 159kg – on flat ground