



BME

BRIE MEDICAL EQUIPMENT

*Local Family
Operated Business*

PRODUCT SHEET

SKU: C8560-A-FX

Peak Care Folding Over Toilet Aid

PRODUCT OVERVIEW

- Superior Corrosion Resistance
- Height Adjustable
- Easy to Assemble
- Easy to Clean
- Compact For Travel



SPECIFICATIONS

Overall Length	529 mm	Seat Height	460 - 620 mm
Overall Height	677 - 827 mm	Unit Weight	3.85 kg
Overall Width	560 mm	Max. Weight Capacity	100 kg
Seat Width	480 mm		