

Shoprider GK10 Crossover

SKU: SC-SH-GK 10

\$3,590



Specifications:

Overall: 1190mm (L) x 585mm (W) x 960mm (H)

Turning Circle: 1330mm

Speed: 9km/h

Maximum range: Up to 27km

Safe climb angle: 9°

Maximum climb angle: 9°

Motor: 1.2Hp high torque motor

Batteries: 2 x 36 Amp minimum (supplied)

Charger: 5 Amp

Tyres (Pneumatic): Ø260mm - Front + Rear



FEATURES :

- Scooter pulls apart for storage and transport
- Comfortable mid back contoured seat-Adjustable steering tiller
- Loop handle with finger operated controls- Rear wheel drive
- Electromagnetic brake
- Battery power indicator
- Split battery box
- Headlights
- Removable rotating / height adjustable seat
- Lift up & adjustable width arm rests
- Lifting handles makes picking up the scooter easy
- Removable shopping basket

Maximum User Weight: 136kg - on flat ground

Product Weight: 78kg - with batteries

Battery Weight: 11kg each

Seat Weight: 16kg

Rear Transaxle Weight: 17kg

Main Chassis Weight: 23kg



