

# Oscar Barwon Dual Motor Lift Chair

SKU: 31056B-2M

**\$2,970**



## ACTIONS



### Single Motor (Wall-saver, Lift Assist)

The single motor wall saver operation of the Barwon allows the chair to move forward while the backrest reclines. This Oscar feature reduces the space that the chair occupies in the room because the reclined backrest does not hit the wall. The Lift, recline and footrest operate from one motor.



### Dual Motor (Independent Recline)

Dual motor option lets the user recline the back without effecting the footrest, allowing a greater range of comfort.

## Features and Benefits:

The Barwon is a Lift-Recline Chair, the first choice for value in a high-quality Australian made lift chair.

Choose from four standard chair size and an extensive range of feature adjustments to achieve a completely personalised fit and ensure comfortable usage.

Every chair is fitted with carefully tested and selected electrical components and can be upholstered in your choice of premium fabric from Oscar's very own Versotela range.

- 130kg weight capacity
- Single motor wall saver
- Controller right side
- 9 volt battery
- 50mm casters
- Comes standard with Velcro attached back cushions
- Emergency battery backup
- Made to measure
- Available in wide range of Versotela fabrics
- 2 Year Warranty on mechanism, electrics and fabric
- 5 Year Warranty on timber frame
- 100% manufactured in Australia by Oscar Furniture



# SPECIFICATIONS



|                  | Seat Height | Seat Width | Seat Depth | Arm Height | Backrest Height | Scapular Height       | Footrest Length |
|------------------|-------------|------------|------------|------------|-----------------|-----------------------|-----------------|
| <b>MINI SIZE</b> |             |            |            |            |                 |                       |                 |
| MIN ADJUSTMENTS  | -30mm       | -150mm     | -150mm     | -60mm      | -180mm          | Linked to back height | -50mm           |
| MEASUREMENT      | 456mm       | 540mm      | 455mm      | 199mm      | 636mm           | 429mm                 | N/A             |
| MAX ADJUSTMENTS  | +20mm       | +50mm      | +75mm      | +60mm      | +150mm          | Linked to back height | +50mm           |
| <b>SIZE A</b>    |             |            |            |            |                 |                       |                 |
| MIN ADJUSTMENTS  | -30mm       | -150mm     | -150mm     | -60mm      | -180mm          | Linked to back height | -50mm           |
| MEASUREMENT      | 476mm       | 540mm      | 490mm      | 180mm      | 647mm           | 440mm                 | N/A             |
| MAX ADJUSTMENTS  | +20mm       | +50mm      | +75mm      | +60mm      | +150mm          | Linked to back height | +50mm           |
| <b>SIZE B</b>    |             |            |            |            |                 |                       |                 |
| MIN ADJUSTMENTS  | -30mm       | -150mm     | -150mm     | -60mm      | -180mm          | Linked to back height | -50mm           |
| MEASUREMENT      | 496mm       | 540mm      | 520mm      | 190mm      | 688mm           | 481mm                 | N/A             |
| MAX ADJUSTMENTS  | +25mm       | +50mm      | +50mm      | +60mm      | +150mm          | Linked to back height | +50mm           |
| <b>SIZE C</b>    |             |            |            |            |                 |                       |                 |
| MIN ADJUSTMENTS  | -30mm       | -150mm     | -150mm     | -60mm      | -180mm          | Linked to back height | -100mm          |
| MEASUREMENT      | 546mm       | 540mm      | 545mm      | 190mm      | 718mm           | 511mm                 | N/A             |
| MAX ADJUSTMENTS  | +50mm       | +50mm      | +25mm      | +60mm      | +120mm          | Linked to back height | +50mm           |