

Ellipse Lite Aluminium Rollator - Medium

SKU: 9346376000512 (Black)

SKU: 9346376000543 (Red)















\$420



Features:

- ▶ AU STD. Tested 492485
- ▶ Adjust Height Tool-Free
- ▶ Adjust Brake Tool-Free
- ▶ Foldable Waterproof Seat
- ▶ Easily Installable Nylon Bag

Length (cm)	 59.5	Folded Width (cm)	 26	Front Castor Diameter (cm)	 24.2/9.5"
Width (cm)	 49	Seat Depth (cm)	 22	Rear Castor Diameter (cm)	 20.3/8"
Max Height (cm)	 97.8	Seat Width (cm)	 45.4	Max Load (kg)	 136
Minimum Height (cm)	 84	Seat Height from the Ground (cm)	 52.3	Net Weight (kg)	 8.2