

# Redgum Compact Side Folding Seat Walker

SKU: RG4401

**\$399**



## FEATURES

- Side folding + collapses for storage / transport
- Strong lightweight construction - Aluminium frame
- User friendly loop brakes with soft brake lock
- Height adjustable handles
- Ergonomic handgrips to reduce pressure
- Comfortable padded seat
- Padded backrest
- Removable carry bag
- Large diameter wheels for excellent mobility
- Flat free rubber tyres
- Swivelling front wheels

## SPECIFICATIONS (All dimensions shown are nominal)

Overall Dimensions (At Lowest Seating): 690mm (L) x 800mm (H) 640mm (Width- Unfolded), 290mm (Folded)

Collapsed: 380mm (L) x 270mm (W) x 620mm (H)

Handle Height (From Floor to top of Handle): 800 - 970mm (8 settings)

Seat: 445mm (W - Between Handles) x 215mm (D) 560mm (From floor to Top of Seat)

Wheels: 8" / Ø200mm (Front - Swivel / Rear - Fixed)

Product Weight: 8.5kg

Max. User Weight: 136kg