

Redgum Bariatric Commode Chair

SKU: B588

\$289



Features and Benefits:

- Allows for toileting away from the toilet for when users cannot make it to the bathroom.
- 280kg user weight capacity
- Plastic seat / backrest / armrests
- Height adjustable
- Removable bowl with carry handle and lid
- Removable splashguard
- Removable backrest
- Nonslip rubber tips
- Lightweight: easy to move
- Easy to clean surfaces
- Bowl, lid and splashguard

Specifications:

Code	B588	Overall Height	700mm (lowest height setting)
Description	Redgum Bariatric Commode Chair	Seat Width	590mm
Overall Width	685mm	Seat Depth	470mm
Overall Depth	620mm	Seat Height	410mm – 575mm
Safe Working Load	140kg	Backrest Pad	400mm (W) x 150mm (H) x 300mm (D)
Product Weight	8.5kg	Maximum User Weight	280kg