

Redgum Basics Commode Chair

SKU: B588

\$235



Features and Benefits:

- Allows for toileting away from the toilet for when users cannot make it to the bathroom.
- Soft padded seat / backrest / armrests
- Removable padded seat overlay
- Easy to clean vinyl upholstery.
- Adjustable seat height
- Flip up seat to access bowl.
- Removable bowl with carry handle and lid
- Nonslip rubber tips
- Lightweight – easy to move.

Specifications:

Code	B580	Overall Height	785mm (lowest height setting)
Description	Redgum Basics Commode Chair with Padded Seat	Seat Width	430mm
Overall Width	610mm	Seat Depth	450mm
Overall Depth	505mm	Seat Height	410mm – 560mm (7 height setting at 25mm increments)
Safe Working Load	140kg	Backrest Pad	420mm (W) x 240mm (H) x 380mm (D) (from seat to top of backrest)
Product Weight	8.5kg	Maximum User Weight	159kg
Width between Armrests	510mm	Product Weight	7kg