


Sovereign S7 Ultimate Care

Foam Surface

SKU: S7SL1000PU (LS)

SKU: S7KSL1000PU (LKS)



Features and Benefits:

Game Changing Foam Innovation

Breathability

- Improved Micro-Climate
- Next generation Ventraflow® provides unrestricted Air Movement which takes heat and moisture away from pressure areas with unique Air Channel System

Adjustable Comfort

- Reversible comfort immersion layer;
- Airoform® cooling Gel memory foam
- Improved body temperature regulation Airoform® Response.
- Highly elasticity for reduced effort while repositioning

Ease of Repositioning

- Ultra-high rebound foam technology makes it easier to roll and reposition.
- The S7 is suitable for High Risk clients, with a weight range between 30 – 250kg.

| Code | Size | Dimensions | Price |
|-------------|------------------|-------------------------------------|---------|
| S7SL1000PU | Long Single | 2030mm (L) x 880mm (W) x 190mm (D) | \$2,197 |
| S7KSL1000PU | Long King Single | 2030mm (L) x 1050mm (W) x 190mm (D) | \$2,397 |