

Redgum MOVO Mobility Flip-Up Scooter Canopy

SKU: SE-MC



\$760

+ Installation \$165 (incl GST)



Features

- Shelter from the elements
- Windshield and Sun visor
- Flip up canopy
- Lightweight
- Easy to fit and remove from scooter
- Other accessories can be fitted behind (extra mount socket)

Minimum fitting requirements:

Mobility Scooter must have the following in order to fit a Movo Canopy:

- Socket under the seat that can take a 25 x 25mm square tube
- Overall width of armrests must be less than 630mm
- Distance from bottom of socket to top of backrest (not including a headrest) must be more than 600mm

Specifications

- Overall: 1500mm L x 750mm W
800mm H
- 870mm Length when folder

