

Shoprider Rocky 4

SKU: 889XLR (Burgundy Red)

SKU: 889XLB (Blue)

\$6,150



BRIEBIE & CABOOLTURE
MEDICAL EQUIPMENT



HIRE AND SALES

Features and Benefits:

- A Powerful 2hp motor
- Large 75amp Gel Batteries
- Front /Rear Suspension & 225kg load capacity
- The Rocky 4 comes with 330mm wheels for added ground clearance
- Airflow Comfort Seat
- Fully rotating lockable seat
- Adjustable steering tiller and multi position lift up armrests.
- High Back Air Flow Comfort Seat
- Loop Steering Tiller
- Adjustable Steering Tiller
- Sliding / Rotating Seat
- Variable Seat Height
- Double Headlights
- Rear Lights & Indicators
- Hazard Lights
- Front & Rear Bump Pads
- Shopping Basket
- Front and Rear Suspension
- High and Low Speed
- Width Adjustable and Lift Up Arm Rests
- L.E.D Battery Level Warning
- Large Battery Capacity

Specifications:

Overall Dimensions: 1430mm (L) x 680mm (W) x 1300mm (H)

Turning Circle: 1800mm

Max. Speed: 10km/h

Maximum range: Up to 40km

Safe climb angle: 9°

Maximum climb angle: 12°

Motor: 2.0Hp high torque motor (4 Pole)

Batteries: 2 x 75 Amp minimum (supplied)

Charger: 8 Amp

Tyres (Pneumatic): Ø330mm - Front + Rear

Product Weight: 146kg - with batteries

Maximum User Weight: 225kg - on flat ground

Warranty:

Ten-Year Warranty:

- Structural Frame

18 Month Warranty:

- Drive Train Assembly

12 Month Warranty:

- Other Parts

Exclusions of Warranty:

1. All transportation costs associated with the submission and return of parts for repair or replacement under warranty (as determined by Shoprider Australia) are the responsibility of the original customer.
2. Normal wear items are not covered under this warranty. These include such items as brushes, tyres, tubes, brake pads, upholstery, seating, and arm rests.
3. Damage caused by:
 - a. Excessive, abusive, misuse, negligence, or accident.
 - b. Improper maintenance, maintenance carried out by non-approved "WARRANTY DEALER".
 - c. Improper operation or storage.
 - d. Any Commercial use.
 - e. Fitment of NON-Manufacturer approved accessories.
Resulting damage attributed to the NON-Manufacturer approved accessory or correction of the unit to operate as intended once the accessory has been removed will not be covered under warranty.
4. Where an authorised dealer has not lodged all the required documentation at the appropriate time, the warranty is void.
5. All warranties provided do not extend to purchasers of used or second-hand equipment.

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement part for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and failure does not amount to major failure.