

Aspire Push Down Seat Walker

SKU: WAF705600RE (Red)

SKU: WAF705600BL (Blue)

\$197



Features and Benefits

- Weight activated brakes ideal for arthritic hands
- Handle height adjusts quickly with hand wheels for optimal sizing
- Rustproof, robust aluminium frame provides durability and stability and whilst remaining very light to manoeuvre and transport
- Easy folding mechanism facilitates storage and transportation
- Anatomically designed, latex-free handles allow comfortable grip especially for users with arthritis
- Curved and padded backrest offers comfortable seated support and flips up/removes easily for folding/storage
- Padded seat flips up easily for access to storage bag (supplied) or for folding
- Swiveling front castor wheels increase maneuverability
- 7 inch diameter solid, puncture-proof wheels provide increased maneuverability over various surfaces
- Available in blue or red
- Easy wipe-clean surfaces

Specifications:

Code	WAF705600RE (red) WAF705600BL (blue)	Seat Width	355mm
Description	Aspire Push Down Seat Walker	Seat Height	580mm
Overall Height	851mm – 962mm	Backrest Height	815mm
Overall Length	652mm	Folded Width	255mm
Overall Width	616mm	Castor Size	7"
Seat Depth	305mm	Unit Weight	7kg