



BME

BRIBIE MEDICAL EQUIPMENT

Local Family

Operated Business

PRODUCT SHEET

SKU: 12325SP

KCare Economy Shower Stool with Plastic Seat

PRODUCT OVERVIEW

- Provides stability and support
- Easy-clean plastic seat
- Corrosion-resistant construction
- Anti-slip feet for safety
- Easy to clean and maintain



SPECIFICATIONS

Overall Height	680 - 810 mm	Seat Height (Adj.)	460 - 600 mm
Seat Depth	300 mm	Seat Width	410 mm
Max. Weight Capacity	110 kg	Material	Powder-Coated Steel Frame