

Aspire Flex Adjustable Seat Walker – 6inch Wheels

SKU: WAF709660CH (Champagne)

\$217



Features and Benefits:

- Super lightweight making it easy to manoeuvre and easy to transport.
- Available in 6" (Champagne) and 8" (Black) wheel variations.
- Robust aluminium frame for improved durability and stability.
- Height adjustable seat height to suit most users.
- Soft padded seat for improved seating comfort.
- Adjustable and lockable brake handles.

Specifications:

Code	SKU: WAF709660CH	Colour	Champagne
Total Length	730mm	Seat Width	360mm
Total Width	570mm	Seat to Floor Height	430-520mm
Total Height	730-830mm	Seat to Footrest	430-520mm
SWL	136kg	Unit Weight	6.3kg