

Aspire Flex Adjustable Seat Walker – 8” Wheels

SKU: WAF709670BK (Black)

\$217



Features and Benefits:

- Super lightweight making it easy to manoeuvre and easy to transport.
- Available in 6" (Champagne) and 8" (Black) wheel variations.
- Robust aluminium frame for improved durability and stability.
- Height adjustable seat height to suit most users.
- Soft padded seat for improved seating comfort.
- Adjustable and lockable brake handles.

Specifications:

Code	WAF709670BK	Colour	Black
Total Length	820mm	Seat Width	360mm
Total Width	585mm	Seat to Floor Height	475-570mm
Total Height	780-970mm	Seat to Footrest	475-570mm
SWL	136kg	Unit Weight	7.3kg