

Aspire Deluxe HD 4 Wheel Scooter – HS898

SKU : SCT671290

\$6,450



Features

The HS898 allows for mobility over challenging terrain. Heavy Duty weight capacity of 181 kg, along with a comfortable captain's seat with angle adjustable back support, depth and swivel adjustment on the seat base and angle adjustable tiller.

A deluxe heavy duty mobility scooter with hill climbing performance!

- Heavy Duty Hill Climber
- Stylish mag wheels
- LED headlights, brake lights and indicators
- Front and rear wheel suspension
- Rearview mirrors
- Adjustable swing away armrests
- Electro mechanical brakes
- Seat Width 520 mm

Range	40km	Safe Working Load	181kg
Battery	2 x 12V 75Ah	Maximum Speed	10km/hr
Total Dimensions (L x W x H)	1440 x 720 x 1360mm	Ground Clearance	100mm
Unit Weight	143kg incl battery	Colour	Black, Blue