

# Merits Pioneer 12

SKU: S8406B-1

# \$5,995



## Colour Options

- ✓ Red
- ✓ Silver
- ✓ Dark Blue
- ✓ Pearl Gold
- ✓ Turquoise



### Features and Benefits:

- Adjustable prestige car like suspension which supplies a smooth ride over varied terrain
- Built-in speed reduction switch which is activated when turning offering exceptional safety
- New state of the art anti-glare LCD screen (dash) with introduction of speedo, odometer, clock, temperature gauge, easy to read more accurate bar style battery gauge. The screen also shows when the scooter is in the half speed mode.
- Dual action parallel tiller adjustment allows the gauges & switches to be visible at all times
- New look dash layout with easy to see & use switches for indicators, head light, horn, hazards & half speed.
- Heavy duty 4 pole motor & larger electronic controller that will deliver all the power you need & provide worry free transport
- Comfortable leather feel contoured high back seat with slide, swivel, recline & armrest height & width adjustment
- New safety reflectors on each side of the scooter for improved visibility
- Power saving LED tail lights, brake lights & indicators at the rear
- Front basket & mirrors standard
- Unique 12 volt USB fitting (for phone charging)

### Warranty Information:

7 years warranty on the frame, 24 months warranty on the drive line, 12 months warranty on other parts & electronic components. (to be checked with Merits Australia)


# Merits Pioneer 12



## Specifications

Speed up to	10kph
Range up to	40-50km (dependent on battery size & terrain)
Turn Radius	161cm
Ground Clearance	13cm
Motor	DC 24V, 400W (Rating)/ 1750W (Max) 4 pole
Controller	Dynamic Rhino 2, 160 amp
Battery	12V/60-80a/h
Charger	8amp off-board
Gradient	12" degrees
Front Wheel	12" Pneumatic
Rear Wheel	12" Pneumatic
Brake	Intelligent, regenerative, electromagnetic
Weight capacity	182kg

## Dimensions



**New Anti Glare LCD Screen**



**Handy USB fitting**



**Adjustable Suspension**



**Powerful 4 Pole Motor**