

Shoprider Rocky 8

SKU: 889AER (Burgundy Red)

SKU: 889AEB (Blue)

SKU: 889AEG (Grey)

\$6,600




Features and Benefits:

- Powerful 4 Pole Motor
- Large High Back Contoured Captain's Seat with Headrest
- Gas Adjustable Steering Tiller
- Lift Up & Adjustable Width Arm Rests
- Sliding Rotating Height Adjustable Seat
- Loop Handle with Finger Operated Controls
- Electromagnetic Brake
- Removable Shopping Basket
- Large Battery Capacity
- Horn
- Deluxe LCD Display Dashboard
- Large Spoke Mag Rims
- Front & Rear Adjustable Suspension
- Front & Rear Lights
- Hazard Lights / Indicators
- High / Low Speed Settings
- Exceptional Ground Clearance

Specifications:

SPECIFICATIONS	SHOPRIDER ROCKY 8
Length	159cm
Width	73cm
Height	128.5cm
Safe Climbing Angle ¹	12 Degrees
Max Climbing Angle ¹	12 Degrees
Maximum Weight Capacity ¹	225kg
Bariatric Suitable	Yes
Charger	8Amp
Motor Size	2.0Hp (4 pole)
Maximum Speed	10km/h
Maximum Range ²	Up to 54km
Wheel Size	33cm
Turning Radius	186cm
Scooter Weight (with batteries)	147kg



Warranty:

Ten-Year Warranty:

- Structural Frame

18 Month Warranty:

- Drive Train Assembly

12 Month Warranty:

- Other Parts

Exclusions of Warranty:

1. All transportation costs associated with the submission and return of parts for repair or replacement under warranty (as determined by Shoprider Australia) are the responsibility of the original customer.
2. Normal wear items are not covered under this warranty. These include such items as brushes, tyres, tubes, brake pads, upholstery, seating, and arm rests.
3. Damage caused by:
 - a. Excessive, abusive, misuse, negligence, or accident.
 - b. Improper maintenance, maintenance carried out by non-approved "WARRANTY DEALER".
 - c. Improper operation or storage.
 - d. Any Commercial use.
 - e. Fitment of NON-Manufacturer approved accessories.
Resulting damage attributed to the NON-Manufacturer approved accessory or correction of the unit to operate as intended once the accessory has been removed will not be covered under warranty.
4. Where an authorised dealer has not lodged all the required documentation at the appropriate time, the warranty is void.
5. All warranties provided do not extend to purchasers of used or second-hand equipment.

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement part for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and failure does not amount to major failure.