

# Aspire Evoke 2 Self-Propelled Wheelchair – 50cm

SKU: MWS449530

**\$997**



## Features and Benefits:

- A lightweight, highly adjustable manual wheelchair
- The contoured padded overlay offers the user increased comfort and is ideal for those who are in their wheelchair for longer periods of time. The mild contouring assists in positioning and pressure distribution, while the fabric allows air movement, increasing breathability. It is easily removed to allow washing when required
- Swing-in / out, height adjustable and removable foot rests
- Strong, Nylon upholstery with a comfortable cushion and padded backrest.
- Depth adjustable seating with rear angle adjustability
- Quick release 24" solid rear tyres

## Specifications:

<b>Seat Depth</b>	400mm - 475mm	<b>Seat Width</b>	500mm
<b>Backrest Height</b>	430mm – 490mm	<b>Seat to Floor Height</b>	440mm – 480mm
<b>Colour</b>	Silver	<b>Material</b>	Aluminium
<b>Safe Working Load</b>	140kg	<b>Unit Weight</b>	15.1kg