

# Freerider Luggie Elite Plus

SKU: FRLUGEP-BK (Black)

SKU: FRLUGEP-NR (Blue)

SKU: FRLUGEP-NR (Red)

# \$6,150



## Features and Benefits:

The Luggie Elite Plus mobility scooter is the same scooter as the Luggie Elite, but with suspension on the front wheels, giving the scooter a softer ride over more bumpy surfaces. The Luggie Elite Plus has all the features of the Elite, including extendable handlebars, bigger Li-ion battery, and armrests.

- Foldable design
- Extendable handlebars
- Armrests as standard
- More powerful motor
- Bigger 10.5 a/h Li-ion battery
- Front-wheel suspension
- Solid tyres
- Narrow body and narrow track
- 145kg weight capacity
- 4.5cm Floor Clearance

## Specifications:

Max. Speed: 6km/h

Maximum range: 15-25km

Batteries: 10.5 a/h Li-ion battery

Charger: 5 Amp

Product Weight: 29kg - with batteries

Maximum User Weight: 145kg

Overall Depth: 99cm

Overall Width: 45cm

Seat Depth: 38cm

Seat Width: 40cm



## **Warranty:**

### **WARRANTY TERMS**

Please keep a note of your serial number.

(This is located on the front frame, under the front chassis cover).

The frame of your scooter is warranted for 36 months, 12 months from the date of purchase against faults arising due to defects in manufacture or materials on electronic components. This warranty does not detract from but is in addition to your legal rights. Parts replaced or repaired under the terms of this warranty will be covered for the balance of the period. This warranty applies ONLY to parts supplied or approved by your Authorized Luggie Dealer. Refer to your warranty sheet for additional details. Due to the individual prescription nature of this product, this warranty is not transferable.

NOTE: Optional extended warranty insurance may be available from your local dealer.

**EXCLUSIONS:** This warranty does not extend to consumable items which may need replacing due to normal wear and tear, namely tires, punctures, lubrication, mats, oil seals, gaskets, screws, upholstery, motor brushes, fuses, batteries, handgrips etc. or to damage to the product caused by misuse, neglect or accident. This warranty does not apply if your Luggie shows signs of excessive wear and abuse or has been modified without the authority of your local dealer.

The following items are not covered by your warranty:

1. Any noise or vibration, which does not affect the quality and function of the scooter.
2. Damage due to lack of maintenance or improper use or operation or storage.
3. Fees for the regular service inspection and maintenance.
4. Any fees incurred as a result of warranty repair, losses or compensation because of the inability to use the scooter (telephone use, shipping, car rental, travel fees etc.)