

# Freerider Luggie Super Deluxe

SKU: FRLUGSD-OB (Blue)

SKU: FRLUGSD-LG (Green)

**\$7,190**



## Features and Benefits:

The Luggie Super Deluxe is the wide track version of the Super Plus, complete with front suspension and solid tyres for those that want more stability and comfort when tackling rugged terrain. The Luggie Super Deluxe comes with a larger seat, and standard 10.5 A/H Li-ion batteries, which can be upgraded to 16.5 A/H Li-ions giving it a much higher range. With the well-designed front suspension and larger tyres on the rear, the ride is smooth and comfortable.

- Compact design
- Powerful motor
- Larger seat
- Wide track stability
- Extendable handles
- Suspension on front axle
- Solid tyres
- 4cm floor clearance
- Max speed 8kph
- 20-25km range
- 10.5 Li-ion batteries (Optional 16.5 A/H)
- 5 amp charger
- Weight capacity 160kg

## Specifications:

Max. Speed: 8km/h

Maximum range: 20-25km

Batteries: 10.5 a/h Li-ion battery

Charger: 5 Amp

Product Weight: 36kg

Maximum User Weight: 160kg

Overall Depth: 99cm

Overall Width: 58cm

Seat Depth: 39cm

Seat Width: 44cm

## **Warranty:**

### **WARRANTY TERMS**

Please keep a note of your serial number.

(This is located on the front frame, under the front chassis cover).

The frame of your scooter is warranted for 36 months, 12 months from the date of purchase against faults arising due to defects in manufacture or materials on electronic components. This warranty does not detract from but is in addition to your legal rights. Parts replaced or repaired under the terms of this warranty will be covered for the balance of the period. This warranty applies ONLY to parts supplied or approved by your Authorized Luggie Dealer. Refer to your warranty sheet for additional details. Due to the individual prescription nature of this product, this warranty is not transferable.

NOTE: Optional extended warranty insurance may be available from your local dealer.

**EXCLUSIONS:** This warranty does not extend to consumable items which may need replacing due to normal wear and tear, namely tires, punctures, lubrication, mats, oil seals, gaskets, screws, upholstery, motor brushes, fuses, batteries, handgrips etc. or to damage to the product caused by misuse, neglect or accident. This warranty does not apply if your Luggie shows signs of excessive wear and abuse or has been modified without the authority of your local dealer.

The following items are not covered by your warranty:

1. Any noise or vibration, which does not affect the quality and function of the scooter.
2. Damage due to lack of maintenance or improper use or operation or storage.
3. Fees for the regular service inspection and maintenance.
4. Any fees incurred as a result of warranty repair, losses or compensation because of the inability to use the scooter (telephone use, shipping, car rental, travel fees etc.)