

Shower Buddy

SKU: SB7E



\$479



Features and Benefits:

Shower Buddy Travel Commode – The premium choice for lightweight, space-saving shower mobility.

It's no secret that a lack of functional space makes bathroom mobility even more challenging. It's common to encounter a compact shower area, particularly when travelling.

The Eco Traveller shower chair from Shower buddy was designed to enable users to shower and use the toilet in bathrooms where space is at a premium

Durable & Reliable: With a lightweight but strong aluminium frame, Eco Traveller can support users up to 120kg, and comes with a 2-year limited warranty.

Portable: Like the name suggests, Eco Traveller is ideal for trips away from home, folding down into a carry bag for ultimate portability.

Functional: Despite its small size, Eco Traveller still enjoys the premium locking caster wheels and complete commode system of our other Showerbuddy products.

Setting up the Eco Traveller is quick and requires no additional tools to assemble:

- Once the unit is assembled, it can be folded down and set up in moments
- Lock the four caster wheels in place and transfer the user onto the chair from the front.
- Release the brake and transfer the user into a roll-in or **accessible shower**.
- After locking the shower chair back in place, the user can then comfortably bathe.

SB7e

ECO TRAVELLER

Folding shower commode

The SB7e travel shower commode chair is a lightweight economical chair designed for accessible or roll-in showers where space is at a premium. It can roll straight over most toilets and into the shower, reducing manual handling and transfers for caregivers.

The SB7e is fully height adjustable, has a soft foam cushion and can also be used as a commode chair, with commode bucket removed from the rear. It has 4 lockable caster wheels and flip up, fold away and height adjustable footrests, has tool-free assembly and quickly folds down to pack away into an accessory carry bag for travel.

Product Features:

- Weight capacity - 120kg
- Rolls over most standard toilets
- Arms - Fixed Height
- Adjustable Footrests - Removable and swing away, height adjustable
- No rust aluminum frame (Grade 6061 T6)
- Stainless steel parts and hardware (Grade 304)
- Caster Wheels - 100mm 4 x locking
- Commode system
- Tool free assembly and easy fold down



Product Dimensions:

- Chair Width - over caster wheels 550mm
- Chair Width - over arms 550mm
- Chair Width ? Inside Arms 450mm
- Chair length - including footrest 840mm
- Backrest Width 448mm
- Backrest Height from Cushion to top 400mm
- Toilet Height Clearance (maximum) 470mm
- Toilet width clearance at narrowest 445mm (at casters)
- Toilet clearance can depend upon the shape of the toilet pan
- Seat Cushion Dimensions 420mm W x 420mm D
- Seat Height from Floor 530mm to 580mm (Adjustable)
- Seat Cushion Opening 175mm W x 335mm L
- Box Dimensions 780mm x 550mm x 210mm
- Box Gross Weight 9.5kg
- Actual Weight of unit 8.9kg